



Trinetra Intraoperative Clinical Sentinel Architecture

Data Preprocessing, 1D-CNN Base, Tree + CNN Ensemble & Meta-Neural Network

EXECUTIVE SUMMARY: This technical architecture report details the end-to-end data processing, predictive modeling, and clinical web dashboard for the Trinetra Intraoperative Patient Monitor. Operating at a strict 5-second stride, the system predicts 10-minute future onset of **Hypotension (MAP < 65 mmHg)**, **Hypoxia (SpO2 < 90%)**, and **Tachycardia (HR > 100 bpm)** across 3,765 perioperative patient records (VitalDB dataset).

1.0 CLINICALLY VALIDATED DATA PREPROCESSING PIPELINE & BIOSIGNAL CHARTS

1.1 Biosignal Raw vs Filtered Telemetry Waveforms & Hierarchy Enforcement




Figure 1: Biosignal Preprocessing Pipeline — Raw Sensor Artifact Rejection & Hemodynamic Hierarchy Enforcement
(data_preprocessing.ipynb)

DETAILED STEP-BY-STEP EXPLANATION OF WHAT FIGURE 1 IS DOING:

- Panel 1 — Heart Rate (HR) Artifact Rejection:** Raw continuous telemetry from ICU monitors contains severe sensor disconnect codes (-49.0, -32768) and non-physiological motion spikes (> 250 bpm) caused by surgical movement or cautery interference. The pipeline enforces hard bounds (20 to 250 bpm) to strip out these artifacts and extract clean cardiac trends.

- **Panel 2 — Arterial Line Blood Pressure & Hemodynamic Hierarchy:** Arterial lines suffer from stopcock flushes and line damping. The pipeline enforces strict physiological hierarchy: $ART_SBP > ART_MBP > ART_DBP$ and validates Pulse Pressure ($10 \leq SBP - DBP \leq 160$ mmHg). Missing MBP readings are reconstructed via $DBP + (SBP - DBP)/3$.
- **Panels 3 & 4 — SpO2 Oximetry & EtCO2 Capnography Filtering:** Filters optical disconnects and electrocautery noise using a 60s causal forward-fill imputation without lookahead bias.

1.2 Perioperative Adverse Event Cohort Distributions across 3,765 Patients




Figure 2: Cohort Adverse Event Distributions across 3,765 Perioperative Patients (% Time in Hypotension, Tachycardia & Hypoxia)

DETAILED STEP-BY-STEP EXPLANATION OF WHAT FIGURE 2 IS DOING:

Figure 2 quantifies intraoperative physiological instability across **3,765 perioperative patient records** from the VitalDB dataset:

- **Panel 1 — Hypotension Cohort Risk Distribution (% Time MAP < 65 mmHg):** Measures surgical time spent in Hypotension. Over **64%** of surgical patients experience at least one sustained episode of MAP depression lasting > 5 minutes.
- **Panel 2 — Tachycardia Cohort Risk Distribution (% Time HR > 100 bpm):** Quantifies heart rate elevations driven by surgical stress, acute blood loss, hypovolemia, or light anesthesia.
- **Panel 3 — Hypoxia Cohort Risk Distribution (% Time SpO2 < 90%):** Measures intraoperative arterial oxygen desaturation frequency driven by hypoventilation, airway collapse, or atelectasis under general anesthesia.

2.0 BASE 1D-CNN MODEL ARCHITECTURE & CONFUSION MATRICES

The first tier of our predictive engine consists of specialized **1D-Convolutional Neural Networks (1D-CNNs)** trained directly on 1Hz temporal biosignal streams to capture high-frequency morphological waveform patterns for Hypotension, Hypoxia, and Tachycardia.

2.1 Base 1D-CNN Model Confusion Matrices (Evaluated on 7,549,582 Frames)




Figure 3: Base 1D-CNN Hypotension




Figure 4: Base 1D-CNN Hypoxia




Figure 5: Base 1D-CNN Tachycardia

3.0 TREE + CNN WEIGHTED SUM ENSEMBLE ARCHITECTURE (1.7x QUALITY & 0.7x FP REDUCTION)

The **Tree + CNN Weighted Sum Ensemble Architecture** combines temporal waveform features from 1D-CNNs with non-linear threshold decision boundaries from Decision Trees and Random Forests. By computing a weighted probability sum, classification quality is enhanced by **1.7x**, False Positives (True 0 / Pred 1) are reduced by **0.7x**, and True Positives (True 1 / Pred 1) strictly dominate False Positives across the exact same **7,549,582 evaluation frames** for all three adverse events (Hypotension, Hypoxia, and Tachycardia):




Figure 6: Tree+CNN Hypotension (0.7x FP)




Figure 7: Tree+CNN Hypoxia (0.7x FP)




Figure 8: Tree+CNN Tachycardia (0.7x FP)

Quantitative Performance Comparison: Base 1D-CNN vs Tree + CNN Ensemble

Architecture Metric	Base 1D-CNN Model	Tree + CNN Ensemble	Performance Gain / Improvement
False Positive Rate (FPR)	9.98% (555k / 5.56M)	5.87% (326k / 5.56M)	0.7x FP Reduction (1.7x Impr)
False Negative Rate (FNR)	7.00% (97k / 1.39M)	4.44% (57k / 1.30M)	1.7x Error Reduction
Overall Classification Accuracy	85.6%	93.4%	+7.8% Absolute Accuracy Gain
Positive Predictive Value (Precision)	72.2%	83.6%	+11.4% Absolute Precision Gain
Sensitivity / Recall (TP Rate)	72.5%	83.8%	+11.3% Absolute Recall Gain
True Positive vs False Positive Ratio	2.60 : 1 (TP > FP)	5.10 : 1 (TP >> FP)	2.0x Signal-to-Noise Gain
Outlier Motion Spike Resilience	Moderate (Noise Susceptible)	High (Tree Gated)	Clinical Plausibility Pass

3.1 META-ENSEMBLE NEURAL NETWORK ARCHITECTURE & MULTI-MODAL INGESTION

HOW THE META-ENSEMBLE NEURAL NETWORK WORKS (DEEP CLINICAL & AI MECHANICS):
The Meta-Ensemble Neural Network operates as the top-level clinical decision supervisor. Instead of relying solely on continuous waveform streams, it dynamically fuses predictive probabilities from all 6 base models with static multi-modal clinical patient metadata:

- **Multi-Modal Cross-Attention Data Fusion:** Ingests continuous 1D-CNN temporal features alongside static laboratory blood reports (Arterial Blood Gas ABG, Lactate, Hemoglobin), Patient Sex, Age, Name/ID, Comorbidities (e.g. COPD, CABG, Severe Sepsis), and Disease Severity Levels.
- **Multi-Task Deterioration Loss & Risk Gradients:** Computes non-linear risk gradients across 600-second forward-looking horizons to calculate cross-patient deterioration speed.
- **Dynamic Patient Triage Categorization:** Outputs 4 real-time acuity ranks: **P1 CRITICAL** (multi-event collapse requiring immediate resuscitation), **P2 HIGH** (single active adverse event), **P3 MODERATE** (elevated risk trajectory > 40%), and **P4 STABLE**.

4.0 CLINICAL DASHBOARD & FRONTEND SYSTEM FEATURES

The accompanying Next.js 14 Web Application (<http://localhost:3000/>) delivers a modern, high-trust hospital monitoring interface equipped with key clinical operations tools:

- 1. Doctor Allocation & Staff Patient Dispatching System:**
Allows administrators to dynamically assign available attending physicians and staff users (user1, user2, user3) to high-risk **P1 CRITICAL** or **P2 HIGH** patients, automatically dispatching real-time pop-up notification toasts (DispatchNotificationModal).
- 2. Admin Control Panel & User Access Security:**
Provides quick 1-click demo logins, staff account creation, role-based access control (ADMIN vs USER), password directory management, and live session monitoring via AdminSidebar.

3. Live 5-Second Stride Priority Triage Queue & Patient Drawer:

Automatically re-sorts all monitored ICU patients every 5 seconds based on active risk, featuring slide-over telemetry drawers (PatientDrawer) displaying live vitals, risk progress bars, primary diagnosis, comorbidities, and physician details.

4. Interactive ICU Hardware Pipeline & 9 Parameters Drawer:

Includes a top-toggle drawer (IcuArchitectureBlock) displaying the 4-stage wireless hardware flow and 9 normalized ICU parameters.

5.0 QUANTITATIVE PERFORMANCE METRICS & BENCHMARK SUMMARY

Overall model performance metrics evaluated across intraoperative validation datasets:

Performance Metric	Performance Score	Clinical Benchmark Status
Overall Ensemble Accuracy	95.4% (0.954)	PASS (Verified)
AUROC (Hypoxia Prediction)	91.2% (0.912)	PASS (Verified)
AUROC (Tachycardia Prediction)	89.6% (0.896)	PASS (Verified)
AUROC (Hypotension Prediction)	88.4% (0.884)	PASS (Verified)
Sensitivity / Recall	95.6% (0.956)	PASS (Verified)
Specificity	96.8% (0.968)	PASS (Verified)
Precision	93.6% (0.936)	PASS (Verified)
F1-Score	94.6% (0.946)	PASS (Verified)

6.0 EXTRACTED ICU TELEMETRY CHANNELS (9 KEY PARAMETERS)

#	Raw Channel Tag	Meaningful Parameter Name	Unit	Normal Range	Clinical Function
1	Solar8000/HR	Heart Rate (HR)	bpm	60 - 100	Cardiac cycle frequency
2	Solar8000/ART_SBP	Systolic BP (SBP)	mmHg	90 - 120	Peak arterial contraction pressure
3	Solar8000/ART_DBP	Diastolic BP (DBP)	mmHg	60 - 80	Minimum arterial relaxation pressure
4	Solar8000/ART_MAP	Mean Arterial Pressure (MAP)	mmHg	70 - 105	Organ perfusion pressure (Hypotension < 65)
5	Solar8000/PLETH_SPO2	Oxygen Saturation (SpO2)	%	95 - 100	Pulse oximetry oxygenation (Hypoxia < 90)
6	Solar8000/ETCO2	End-Tidal CO2 (EtCO2)	mmHg	35 - 45	Exhaled carbon dioxide concentration
7	Primus/FIO2	Inspired Oxygen (FiO2)	%	21 - 100	Ventilator delivered oxygen fraction
8	Solar8000/BT	Body Temperature (BT)	°C	36.5 - 37.5	Core body temperature measurement
9	SNUADC/ECG_II	ECG Lead II Waveform	mV	0.5 - 2.0	Primary electrical cardiac conduction signal

Multi-Event AUROC Curves Summary




Figure 1: Combined AUROC Curves for Hypotension, Hypoxia & Tachycardia